

Proact Metal Max™+ Combi 50

Laryngoskopset - enpatientmodel



- “Alt i et” laryngoskop af enpatientmodel. ”Test-in-bag” forpakning muliggør funktionstest inden anvendelse
- Oplades med 6V batteri
- Superljus, 25 Lumen LED-lyskilde gir minst 10 000 lux
- Opfylder kravet enligt ISO 7376:2020

- Unikt bajonetklipps-batterilås gør det muligt at tage batteriet væk på mindre end 3 sekunder
- Ergonomisk håndtag i blød silikone for et udmærket greb
- Holdbarhed: 3 år fra fremstillingsdato
- Material: rustfrit stål, nylon, aluminium, silikon, TPE (termoplastisk elast/gummi) og polykarbonat

Bestillingsoplysninger

Artikelnummer	Håndtag	Blad design	Størrelse*	Blad længde, mm	Bladets bredde, mm	Pakn
MMCB50MAC0	Paediatric	MAC 0	Neo	61,0	8,5	10 st/pkt
MMCB50MAC1	Paediatric	MAC 1	Small child	75,0	11,5	10 st/pkt
MMCB50MAC2	Adult	MAC 2	Child	93,0	12,5	10 st/pkt
MMCB50MAC3	Adult	MAC 3	Adult	110,0	13,5	10 st/pkt
MMCB50MAC4	Adult	MAC 4	Large adult	135,0	14,0	10 st/pkt
MMCB50MAC5	Adult	MAC 5	X-large adult	156,0	14,0	10 st/pkt
MMCB50MIL00	Paediatric	Miller 00	Premature	44,0	10,5	10 st/pkt
MMCB50MIL0	Paediatric	Miller 0	Neo	55,0	10,5	10 st/pkt
MMCB50MIL1	Paediatric	Miller 1	Small child	81,0	10,5	10 st/pkt
MMCB50MIL2	Adult	Miller 2	Child	131,0	13,0	10 st/pkt
MMCB50MIL3	Adult	Miller 3	Adult	171,0	13,0	10 st/pkt
MMCB50CP00	Paediatric	Cardiff PRO 00	Premature	55,0	8,0	10 st/pkt
MMCB50CP0	Paediatric	Cardiff PRO 0	Neo	71,0	9,0	10 st/pkt
MMCB50CP1	Paediatric	Cardiff PRO 1	Small child	77,0	10,0	10 st/pkt