

Proact Metal Max™+ Combi 90

Laryngoskopset - enpatientmodel



- “Allt i ett” laryngoskopset av enpatientmodel
- Inkluderer 6V batteri med 3 års holdbarhed fra fremstillingsdatoen og lang batterilevetid, 5 timers kontinuerligt LED-lys
- Förpackad så att det går att funktionstesta innan användning, ”Test-in-bag”
- Miljøvenlig alkalisk batteriteknik så batteriet er enkelt at fjernes til genbrug
- Findes tillgængelig i MAC, Miller, Lever Tip og Cardiff PRO utförande
- Material: Rustfrit stål, Aluminium, Polykarbonat

Artikelnummer	Håndtag	Blad design	Størrelse*	Blad længde, mm	Blad bredde, mm	Pakn
MMCB90MAC0	Paediatric	MAC 0	Neo	61,0	8,5	10 st/pkt
MMCB90MAC1	Paediatric	MAC 1	Small child	75,0	11,5	10 st/pkt
MMCB90MAC2	Paediatric	MAC 2	Child	93,0	12,5	10 st/pkt
MMCB90MAC3	Adult	MAC 3	Adult	110,0	13,5	10 st/pkt
MMCB90MAC3S	Stubby	MAC 3	Adult	110,0	13,5	10 st/pkt
MMCB90MAC4	Adult	MAC 4	Large adult	135,0	14,0	10 st/pkt
MMCB90MAC4S	Stubby	MAC 4	Large adult	135,0	14,0	10 st/pkt
MMCB90MAC5	Adult	MAC 5	X-large adult	156,0	14,0	10 st/pkt
MMCB90MIL00	Paediatric	Miller 00	Premature	44,0	10,5	10 st/pkt
MMCB90MIL0	Paediatric	Miller 0	Neo	55,0	10,5	10 st/pkt
MMCB90MIL1	Paediatric	Miller 1	Small child	80,5	10,5	10 st/pkt
MMCB90MIL2	Paediatric	Miller 2	Child	131,0	13,0	10 st/pkt
MMCB90MIL3	Adult	Miller 3	Adult	170,0	13,0	10 st/pkt
MMCB90MIL4	Adult	Miller 4	Large adult	182,0	17,0	10 st/pkt
MMCB90LT2	Paediatric	Lever Tip 2	Child	93,0	12,5	10 st/pkt
MMCB90LT3	Paediatric	Lever Tip 3	Adult	110,0	13,5	10 st/pkt
MMCB90LT4	Adult	Lever Tip 4	Large adult	135,0	14,0	10 st/pkt
MMCB90CP00	Paediatric	Cardiff PRO 00	Premature	55,0	8,0	10 st/pkt
MMCB90CP0	Paediatric	Cardiff PRO 0	Neo	71,0	9,0	10 st/pkt
MMCB90CP1	Paediatric	Cardiff PRO 1	Small child	77,0	10,0	10 st/pkt

* Størrelse enligt ISO 7376 standard